

## PRE-OWNED AUTOREF/KERATOMETERS

Refurbished, like new, patient ready!



### Nidek OPD-Scan III Wavefront Aberrometer

The OPD-Scan III is an autorefractor, keratometer, pupillometer (up to 9.5 mm), corneal topographer, and integrated wavefront aberrometer. The OPD-Scan III completes 20 diagnostic metrics in less than 10 seconds per eye (including angle kappa, HOAs, average pupil power, RMS value, and point spread function). Easy alignment and automatic capture of wavefront aberrometry data ensures accurate readings. Wavefront aberrometry data is gathered from available zones up to a 9.5 mm area, adding the capability to provide for calculation of mesopic refractions. Blue light, 33 ring, placido disc topography is gathered in one second. Mapping methods include OPD, visual acuity corneal topography/topographer, and more.

**SAVE \$2,000**

One-year warranty and free shipping use code MYFREE



### Nidek Tonoref II ARK/Tono

- Combines the technologies of autorefractor, keratometer, and non-contact tonometer
- Smaller footprint, to save even more space,
- Uses SLD (Super Luminescent Diode) and a highly sensitive CCD device to provide Zonal Ring-Image Technology.
- Accurately measures patients with cataracts, corneal opacities, IOLs, and post LASIK.

**SAVE \$2,000**

One-year warranty and free shipping use code MYFREE



Topcon KR 8000 PA The three in one

- **It is an Autorefractor, Keratometer, and Topographer.**

- Pupil dilation to 2.0 mm.
- 3D auto-alignment system for easy operation.
- Completes three readings of each eye.
- Displays corneal mapping very quickly.

The wider area of 1.5 mm ~ 9.2 mm can be measured using 10 Placido rings.

**SAVE \$1,000**

One-year warranty and free shipping use code MYFREE



### Topcon KR 8900

- Rotary Prism measurement for enhanced accuracy
- Refractive and keratometric measurements through pupils as small as 2 mm.
- An easy-to-read color LCD display
- Autofire mode for fast, reliable data collection
- Built-in printer

**SAVE \$1,000**

One-year warranty and free shipping use code MYFREE



### Topcon KR 8800

- Color LCD monitor simultaneously displays the readings of both right and left eyes.
- Allows measurements to be obtained on patients with up to 10D of astigmatism.
- Rotary Prism Technology for unmatched accuracy and reliability.
- The scenic color chart for easy eye fixation.

- Completes three readings of each eye.
- Measures the corneal diameter either during testing of the patient or later by recall.

**SAVE \$1,000**

One-year warranty and free shipping use code MYFREE





Topcon KR 8000

- Rotary Prism Technology for unmatched accuracy and reliability.
- Capable of obtaining perfectly accurate keratometry and/or refraction measurements pupils as 2.0 mm.

**SAVE \$1,000**

One-year warranty and free shipping use code MYFREE



Marco-Nidek ARK 760A

- With both visual acuity and subjective refinement.
- Electronic elevation, automatic alignment, focusing, printing, high-speed measurement.

One-year warranty and free shipping use code MYFREE

**For a Bigger Return on Your Investment**

**Call, text, WhatsApp NOW: 949-362-1150**

[info@usedeyeequipment.com](mailto:info@usedeyeequipment.com)