

BEST DEALS

PRE-OWNED CORNEAL TOPOGRAPHERS


Refurbished, like new, patient ready!



Nidek OPD-Scan III Wavefront Aberrometer

The OPD-Scan III is an autorefractor, keratometer, pupillometer (up to 9.5 mm), corneal topographer, and integrated wavefront aberrometer. The OPD-Scan III completes 20 diagnostic metrics in less than 10 seconds per eye (including angle kappa, HOAs, average pupil power, RMS value, and point spread function). Easy alignment and automatic capture of wavefront aberrometry data ensures accurate readings. Wavefront aberrometry data is gathered from available zones up to a 9.5 mm area, adding the capability to provide for calculation of mesopic refractions. Blue light, 33 ring, placido disc topography is gathered in one second. Mapping methods include OPD, visual acuity corneal topography/topographer, and more.

One-year warranty and free shipping use code MYFREE



Humphrey Zeiss Atlas

The Atlas System Model 995 and 900 features include advanced Arc-Step Algorithm, true elevation data, and aspheric surface reconstruction.



Oculus Pentacam

Equipped with an automatically rotating Scheimpflug camera the Pentacam performs a complete measurement of the anterior eye segment in less than 2 seconds.

OUT OF STOCK

For a Bigger Return on Your Investment

Call, text, WhatsApp: 949-362-1150

info@usedeyeequipment.com

