



Right Medical Speedy K

Reconditioned, Calibrated, Patient Ready



Speedy K by Right Medical, offers advanced engineering, rapid acquisition of data, and the medical retro illumination mode to make the Speedy K autorefractor the most advanced autorefractor available.

ACCURATE: Advanced engineering, rapid acquisition of data, and the medical retro illumination

mode makes the Speedy autorefractors the most advanced autorefractor available.

EASY: Large, 5.5 inch high quality monitor allows for easy alignment and focusing.

FAST: Readings are taken the instant that alignment and focus is achieved, and terminated when accurate data is gathered. Measurements are taken at rapid speed, making the instrument perfect for patients who require fast measurement acquisition.

Simultaneous refraction and keratometry measurement can be taken in 0.35 seconds.

**For a Bigger Return on Your Investment**

**Call, text, WhatsApp: 949-362-1150**

**[info@usedeyeequipment.com](mailto:info@usedeyeequipment.com)**